

Workshop and Training Program on Sampling and Detection Methods Applied to Transgenic Crops

November 17 – 19, 2011, NIN, Hyderabad, India

Quantitative PCR for GMO Detection


E. Pearce Smith
Eurofins|GeneScan

NIN, Hyderabad
November 17-19, 2011

1

DNA-based Analytical Process

2




In practice:

- Efficiency is not 100%
- Depletion of PCR reagents

Therefore:

- Exponential increase is limited
- Linear increase follows exponential and plateau occurs

4



What is Real-Time PCR?

A PCR system in which we can monitor the amplification reaction as it is occurring

Real-Time PCR incorporates the ability to directly measure and quantify the reaction while amplification is taking place

6

Workshop and Training Program on Sampling and Detection Methods Applied to Transgenic Crops

November 17 – 19, 2011, NIN, Hyderabad, India

Which detection strategies are available ?

Intercalating dyes

Hybridisation probes

7

Intercalating Dyes - SYBR green

Binds to double-stranded DNA (dsDNA).
Fluoresces only when bound to dsDNA.
As PCR progresses [dsDNA] increases, fluorescent signal in sample increases.

Repeat

8

Hybridization Probes

TaqMan Fluorogenic Probe

Reporter (Fluorescein) Quencher (Rhodamine)

Energy Transfer

5' 3'

Fam R Vic Phosphate Group

Laser Excitation

9

TaqMan PCR Chemistry

R = Reporter
Q = Quencher

10

TaqMan Minor Groove Binder (MGB) probes

AGGCCTTGAGAGATAT NFQ

R MGB

AGGCCTTG Q GATAT

GCTACACAGTCCGGAACCTCTATAGCATCAGAC

- MGB are peptides derived from natural antibiotics that bind DNA in the minor groove stabilizing duplex DNA/Probe
- No Fluorescent Quenchers (NFQ) lose excitation energy by emission of heat (instead of light as in the case of TAMRA).
- Advantages: increased stability and specificity, lower background

11




Which is the best phase to extract quantitative information?

96 PCR replicates of identical samples have very different individual efficiencies at the end of the reaction

12


Workshop and Training Program on Sampling and Detection Methods Applied to Transgenic Crops

November 17 – 19, 2011, NIN, Hyderabad, India



GMO quantitation by real-time PCR: what do we need ?

- Standard Curve prepared with material at known GM content, commonly Certified Reference Materials
- Amplification of two targets:**
 - Transgene
 - Endogenous species-specific gene



The "Standard Curve" method

Two separate Standard Curves: **transgene** and **endogenous gene**

GMO quantification

GM TARGET	Copy Number
Standard 1	2000
Standard 2	1000
Standard 3	500
Standard 4	250
Standard 5	125
Unknown 1	1853
Unknown 2	654
Unknown 3	312

Endogenous gene quantification

ENDOGENOUS	Copy Number
Standard 1	100000
Standard 2	50000
Standard 3	25000
Standard 4	12500
Standard 5	6250
Unknown 1	87500
Unknown 2	73280
Unknown 3	93989

Workshop and Training Program on Sampling and Detection Methods Applied to Transgenic Crops

November 17 – 19, 2011, NIN, Hyderabad, India

Percentage of Transgenic Material

The percentage of transgenic material (DNA) is determined as being the **ratio** of transgenic to total DNA quantities

$$\text{GM DNA \%} = \frac{\text{DNA Q transgenic system}}{\text{DNA Q endogenous system}} \times 100$$

Sample	Transgene Copy Number	Endogenous copy number	% GMO
Unknown 1	1853	87500	2.12
Unknown 2	654	73280	0.89
Unknown 3	312	93989	0.33

Plate layout for GMO quantitation: example of "simplex" assay

S1	S1	S1	S2	S2	S2	S3	S3	S3	S4	S4	S4	G M O
A	A	A	A	A	A	B	B	B	B	B	B	
C+	C+	C+	C-	C-	C-	NTC	NTC	NTC				
S1	S1	S1	S2	S2	S2	S3	S3	S3	S4	S4	S4	E N D O G E N O U S
A	A	A	A	A	A	B	B	B	B	B	B	
C+	C+	C+	C1	C1	C1	NTC	NTC	NTC				

C+ Positive DNA target control

C- Negative DNA target control

NTC Amplification reagent control

What and how do we measure ?
(Some) Critical factors affecting quantitative measurements in GMO analysis.

Target copy number → Discrepancies in GM target copy number between standard and unknown samples may lead to at least 100% over/underestimation of the GM content

Stability of reference gene → Reference gene should be tested for stability of amplifiable target copy number across species (i.e. allelic variation)

PCR efficiency/inhibition → Sub-optimal PCR efficiency due to inhibition may greatly affect linearity of the system and thus correctness of measurements.

Main Critical performance indicators for quantitative PCR methods in GMO analysis

Specificity >>> Primer design and testing, choice of unique target sequences. Occurrence of false positive results

Sensitivity >>> PCR design and optimisation, quality of DNA template (degree of DNA degradation, presence of inhibitors, etc.). Limit of Detection and false negative results.

Precision and trueness >>> Method optimisation, inter-laboratory transferability. Limit of Quantification

What does '% GMO' mean ?
What is the Unit of Measurement ?

weight : weight ?
seed : seed ?
mass : mass ?
protein : protein ?
DNA : DNA ?
..... ??

'GMO' Quantification by PCR - A Ratio, not an Absolute Amount

	Corn DNA*	'GMO' DNA*	Result on test report
Corn	80,000	800	1 %
Corn gluten	8,000	80	1 %
Corn starch	800	8 ?	Not quantifiable
Maltodextrin	80	0.8 ??	Not quantifiable
Glucose syrup	0 ?	???	Not quantifiable

*Approximate relative numbers of target DNA molecules that can typically be extracted from a sample and analyzed in quantitative PCR reactions

- Percentages represent the *relative* composition of the DNA solution extracted from a sample; amount of biotech DNA relative to the amount of the plant species DNA.
- This percentages is not a measure of the absolute amount of biotech DNA in a sample.


ILSI
International Food Biotechnology Committee

Workshop and Training Program on Sampling and Detection Methods Applied to Transgenic Crops


November 17 – 19, 2011, NIN, Hyderabad, India

Testing for % weight / weight?


- Stacked events result in an analytical overestimate of the percentage of biotech corn as compared to the actual w/w percentage of biotech corn in the sample.




1 kernel in 100
PCR result = 1 %
w/w = 1 %




1 kernel in 100
PCR result = 1 %
w/w = 1 %



1 kernel in 100
PCR result = 2 %
w/w = 1 %

 ILSI
International Food Biotechnology Committee
11/29/2011

Adventitious presence of stacked events




Bulk samples

~~Homogenization~~


DNA Extraction

PCR tests


Adapted from Hiroshi Akiyama, NIHSS, Japan

 ILSI
International Food Biotechnology Committee

Analysing for combined events




- No PCR or other molecular approach can distinguish between the presence of a low percentage of stacked events in bulk grain and the presence of a mixture of the two or more individual events that comprise the stack.

 ILSI
International Food Biotechnology Committee


Analysing for combined events

- Analysis of Single seeds or sampling approaches have been proposed:
 - Akiyama et al. proposed testing individual grains
 - difficult and expensive
 - only feasible for high % GM
 - The International Seed Testing Association (ISTA) has proposed a statistical approach for seed lots
 - Only effective for seed, for limited numbers of stacks and at low concentrations

 ILSI
International Food Biotechnology Committee

Challenges & Technical Limitations of PCR

- Is sample representative
- DNA quality and abundance
- Genetic makeup of the sample (i.e. stacked events)

 ILSI
International Food Biotechnology Committee
11/29/2011